## **Redox Titration in an Acidic Medium**

## Laboratory Simulation for the Titration between Potassium Permanganate and Iron (II)

$$5 \text{ Fe}^{2+}(aq) + \text{MnO}_4(aq) + 8 \text{ H}^+(aq) \rightarrow 5 \text{ Fe}^{3+}(aq) + \text{Mn}^{2+}(aq) + 4 \text{ H}_2\text{O}(1)$$

1) Using the reaction above, fill out the table below:

Element	Oxidation Number as a Reactant	Oxidation Number as a Product	Gaining or Losing Electrons?	Oxidized or Reduced?
Fe				
Mn				

2) How many electrons are transferred in this redox titration?

 If the titration is known to be complete when it changes color, is this the end point or equivalence point of the titration? End point or Equivalence Point

- 4) Google "redox titration simulation".
- 5) Click on REDOX Titration in Acidic Solution Computer Simulation
- 6) Click on the link <u>http://pages.uoregon.edu/tgreenbo/redox.html</u>
- 7) You may have to Agree to the disclaimer.
- 8) You may have to enable Adobe Flash Player by clicking the right-hand side of the address puzzle piece, click Manage, unclick "Block Sites from Running Flash (recommended)", and then refresh your page.
- 9) Once the simulation is up, click on 1. Select the Reaction.
- 10) Click on the Oxidizing Agent KMnO<sub>4</sub> and the Reducing Agent Fe<sup>2+</sup>.
- 11) What purple solution was filled in the buret?
- 12) What is the Molarity of KMnO<sub>4</sub>?
- 13) Start your titration by adding the KMnO<sub>4</sub>. Add the KMnO<sub>4</sub> by 2. Push Slider Up to Add a Volume of KMnO<sub>4</sub>. Do not put too much in. Add about 1-3 mL each slide. Keep adding until you can start to see a color change. Once you see a color change, you may want to add using the Dropwise button. Do not go past your end point! Once the color totally changes, you are done your titration.
- 14) What was the Volume of KMnO<sub>4</sub> you added?
- 15) Calculate the number of moles of KMnO<sub>4</sub> added.
- 16) What is the molar ratio of moles of  $MnO_4^-$ : moles of  $Fe^{2+}$ ?
- 17) Calculate the number of moles of Fe<sup>2+</sup> that was reacted.
- 18) What was the Volume of Fe<sup>2+</sup> initially?
- 19) Calculate the Molarity of Fe<sup>2+</sup> (round to 2 significant digits).
- 20) Enter your Molarity of Fe<sup>2+</sup> into 3. After Titration, Calculate and Enter Molarity of Fe<sup>2+</sup>.
- 21) Did you calculate the correct Molarity? Yes or No
- 22) If you went past the end point or if you go the Molarity incorrect, you can always Repeat the Experiment by Pressing Repeat. If you got it correct, click Reset and do the next Oxidizing Agent.

## Laboratory Simulation for the Titration between Potassium Dichromate and Tin (II)

 $3 \operatorname{Sn}^{2+}(aq) + \operatorname{Cr}_2 \operatorname{O}_7^{2-}(aq) + 14 \operatorname{H}^+(aq) \rightarrow 3 \operatorname{Sn}^{4+}(aq) + 2 \operatorname{Cr}^{3+}(aq) + 7 \operatorname{H}_2 O(1)$ 

1) Using the reaction above, fill out the table below:

Element	Oxidation Number as a Reactant	Oxidation Number as a Product	Gaining or Losing Electrons?	Oxidized or Reduced?
Sn				
Cr				

Molarity of K <sub>2</sub> Cr <sub>2</sub> O <sub>7</sub>	
Volume Used of K <sub>2</sub> Cr <sub>2</sub> O <sub>7</sub>	
Moles Used of K <sub>2</sub> Cr <sub>2</sub> O <sub>7</sub>	
Molar Ratio of Moles of $K_2Cr_2O_7$ : Moles of $Sn^{2+}$	
Moles of Sn <sup>2+</sup>	
Volume of Sn <sup>2+</sup>	
Molarity of Sn <sup>2+</sup>	

## Laboratory Simulation for the Titration between Iodine and Thiosulfate Ion

 $I_2(l) + 2 S_2 O_3^{2-}(aq) \rightarrow S_4 O_6^{2-}(aq) + 2 I^-(aq)$ 

1) Using the reaction above, fill out the table below:

Element	Oxidation Number as a Reactant	Oxidation Number as a Product	Gaining or Losing Electrons?	Oxidized or Reduced?
Ι				
S				

Molarity of I<sub>2</sub>

Volume Used of I<sub>2</sub>

Moles Used of I<sub>2</sub>

Molar Ratio of Moles of  $I_2$ : Moles of  $S_2O_3^{2-}$ 

Moles of S<sub>2</sub>O<sub>3</sub><sup>2-</sup>

Volume of S<sub>2</sub>O<sub>3</sub><sup>2-</sup>

Molarity of S<sub>2</sub>O<sub>3</sub><sup>2-</sup>